GRANBERG 0022E

REF 210.0022E

EN: User Manual

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SV: Bruksanvisning

DA: Brugervejledning

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User Manual issue date: 25.11.2025

Head office: GRANBERG AS,
Bjoavegen 1442, 5584 Bjoa, Norway.
Phone: +47 53 77 53 00
E-mail: post@granberg.no

Swedish office: GRANBERG SVERIGE AB,
Schubergsvägen 20, 311 74 Falkenberg, Sweden.
Phone: +46 (0)346 124 25
E-mail: post@granberg-ab.se

GRANBERG®

granberggloves.com

GRANBERG 0022E

REF 210.0022E

Explanation of symbols and pictograms used

	Consult instructions for use		Corrugated cardboard
	Manufacturer		Non-corrugated cardboard
	Date of manufacture		Paper
	Expiry date		Low density polyethylene
	Lot number		Medical Device
	Fragile, handle with care		Unique Device Identifier
	Article number		
	Keep away from sunlight		
	Keep dry		
	Temperature limit		
	Non-sterile		
	Do not re-use		
	Suitable for contact with foodstuffs. Check the Food Declaration of Compliance.		

ENGLISH

210.0022E – GRANBERG SURGICAL DISPOSABLE FACE MASKS TYPE IIR

Description

Disposable surgical face masks Type IIR are made of non-woven polypropylene. Material in 3-layers with elastic loop for ears. The face masks are 100% latex free, >99% bacterial filtration efficiency.

The face masks are made of following material:

Outer layer: Polypropylene.

Second layer: Meltblown fabric.

Inner material: Polypropylene.

Ear loops: Polyurethane

Nose clip material: Polypropylene and steel.

Medical Device Regulation

This product is classified as a Class I medical device and meets the provisions of the Medical Device Regulation (EU) 2017/745. Non-sterile, single use.

Intended purpose: medical face masks are designed to prevent transmission of infective agents, blood and/or bodily fluids from the person wearing them (including but not limited to healthcare personnel and patients) onto others.

Performance

This product meets the requirements of the European Standard EN 14683:2025:

Surgical face masks, type IIR.

Test	Type I *	Type II	Type IIR
Bacterial filtration efficiency (BFE), (%)	≥ 95	≥ 98	≥ 98
Differential pressure (Pa)	≤ 200	≤ 200	≤ 300
Splash resistance pressure (kPa)	Not required	Not required	≥ 16,0
Microbial cleanliness (cfu/g)	≤ 30	≤ 30	≤ 30

* Type I medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations. Type I masks are not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.

Face masks comply with Regulation (EC) No. 1935/2004 and Regulation (EC) No. 2023/2006, and according to Commission Regulation (EU) No. 10/2011, they are suitable for handling all types of foods for short periods of time.

Material Performance

The performance of the face masks quoted above is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, mechanical wear etc.

User guide

1. Wash your hands thoroughly before you put on the mask.
 2. Hold the mask by its edges in front of you, with the white side inwards and the nose clip upwards. Pull a rubber band behind each ear.
 3. Pinch and stroke your fingers along the nose clip a couple of times to form it according to the shape of your nose ridge, preventing air from entering between the mask and your nose.
 4. Pinch the nose clip while pulling the lower part of the face mask below your chin to achieve a better and tighter fit.
 5. Do not touch the mask while wearing it, and only hold it by the rubber bands when removing it. Discard the mask in a proper manner.
- Do not use the mask for more than 4 hours. Remove it if breathing becomes difficult, it becomes damp or the fit of the mask becomes loose.

Storage

This product should be stored in a dry, ventilated environment where the temperature is between -5°C to +30°C. Keep away from direct sunlight, fire and flammable substances. Preferably stored in original packaging.

The shelf life for products stored as recommended is mentioned on each package. Lifetime cannot be specified as it depends on its use. It is the responsibility of the user to determine whether the face mask is suitable for intended use.

Reporting of serious incidents

In case of any serious incident has occurred in relation to the device, please report to the manufacturer and the competent Authority.

Use limitations and Disposal

- The used raw materials are chosen to minimize the risk of allergy and hypersensitivity, but a guarantee for allergic reactions cannot be given. If any allergic reaction occurs, seek medical advice or visit a doctor immediately.
- EU Declaration of Conformity is accessible at www.granberg.no/catalog/210.0022E
- For single use only. Do not re-use.
- Before usage, inspect face mask for any defect or imperfections.
- Dispose according to hospital policy and local regulation.
- The manufacturer cannot be held responsible for improper use of the face mask.

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