Instructions for Use—English

Model 1000USB-C USB Interface Adapter

Indications for Use

Nonin's 1000USB-C interface adapter is designed for connecting Nonin handheld pulse oximeter to a Continua-compatible Application Hosting Device (AHD) or personal computer (PC). Data is transferred between the pulse oximeter and the Continua AHD or PC when you use this product with a Continua Certified™ AHD or PC that supports the ISO/IEE 11073-10404:2010 (Device Specialization-Pulse Oximeter).

Warnings

- Setup of the device should be done in the order listed below to prevent the risk of shock to patient/operator.
- The operator shall not touch the patient and the accessible pins on the accessory simultaneously to prevent the risk of shock to patient/operator.
- Accessory equipment connected to the analog and digital interfaces must be in compliance with the respective nationally harmonized IEC standards (i.e., IEC 60950 for data processing equipment, IEC 61010-1 for laboratory equipment, and IEC 60601-1 for medical equipment).
- To prevent improper performance and/or patient injury, verify compatibility of the monitor, sensor(s), and accessories before use.
- To avoid patient injury, use only Nonin-branded accessories.

Note: Refer to Nonin's handheld pulse oximeter operator's manuals for additional contraindications, warnings, and cautions.

Notes:


- The 1000USB-C interface adapter is designed for use with all Nonin handheld pulse oximeters.
- The 1000USB-C interface adapter may be used for real-time data collection only. The Nonin memory playback feature is not supported at this time.
- The Application Hosting Device (AHD) or PC, AHD application, and corresponding USB Personal Health Device Class (PHDC) host driver are the responsibility of the AHD vendor. Contact the AHD vendor for product support on these items.
- This device does not indicate connection problems associated with the Continua AHD or PC.
- Users and/or patients should report adverse events involving their Nonin device to Nonin Medical, Inc. and the competent authority of the EU Member State in which the user and/or patient is established, if applicable.
- The device's expected service life is 1 year.

Symbols:

Symbol	Definition of Symbol
	Consult Instructions for Use.
	Follow Instructions for Use.
	CE Marking indicating conformance to EC Directive No. 93/42/EEC concerning medical devices.
	Type BF Applied Part
	Indicates separate collection for electrical and electronic equipment (WEEE).
	Protected against vertically falling water drops when enclosure is tilted up to 15 degrees and ingress of solid foreign objects greater than or equal to 2.5 mm (0.1 in.) in diameter per IEC 60529.
	Continua Certified™ signifies that this product has been tested and proven to be interoperable with other products that carry the "Continua Certified" symbol.
	USB Connector Port
	Medical Device
	Importer
	Unique Device Identifier
	Country of Manufacture
	Manufacturer
	Date of Manufacture
	Lot Number
	Serial Number
	Quantity
	Catalogue Number
	Authorized Representative in the European Community
	Handle with care
	Keep dry
	Storage/Shipping Temperature Range

Setup:

- Follow the instructions provided by your Continua AHD vendor to set up your Continua AHD or PC.
- Connect the sensor connector of the pulse oximeter to the "Oximeter" adapter port (1).
- Connect the sensor to the "Sensor" adapter port (2).
- Connect the "→←" adapter port (3) to the USB port on your Continua AHD or PC.
- Follow the instructions provided by your Continua AHD vendor for collecting data from a Continua Certified™ pulse oximeter.

Connecting the Device into a Medical System

Incorporating the device into a medical system requires the integrator to identify, analyze, and evaluate the risks to patient, operators, and third parties. Subsequent changes to the medical system after device integration could introduce new risks and will require additional analysis. Changes to the medical system that must be evaluated include:

- Changing the system configuration
- Adding devices to or disconnecting devices from the system
- Updating or upgrading equipment connected to the system

Issues resulting from user-initiated system changes may include corruption or loss of data.

NOTES:

- When using the serial port to connect the device to other equipment, follow the manufacturer's clearing instructions.
- Verify all equipment connected to the device is suitable for the patient's environment.

Care and Maintenance

Wipe the device with a soft cloth dampened with a mild detergent or 10% bleach solution. Do not use undiluted bleach or any cleaning solution other than those recommended here, as permanent damage could result. Dry with a soft cloth or allow to air dry.

Specifications

Temperature:
Operating: -20 °C to 50 °C (-4 °F to 122 °F)
Storage/Transportation: -40 °C to 70 °C (-40 °F to 158 °F)
Humidity:
Operating: 10% to 90% non-condensing
Storage/Transportation: 10% to 95% non-condensing
This product complies with ISO 10993-1, Evaluation Biological of Medical Devices Part 1: Evaluation and Testing.

Warranty

90 days from the date of delivery.
www.nonin.com/warranty

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